


1kb Ladder

M1181/M1182
50 /50 |5



3-5 ×l/

1kb ladder

8 cm 1%

1l TAE 7 V/cm 45 min

-20

100 ng/5×l

40 ng/5×l

1kb ladder 10 DNA 500 bp 10 kb DNA
2000 bp 5000 bp

104 ng/×l

bp
500 1,000 1,500 2,000 3,000 4,000 5,000 6,000 8,000 10,000

DNA Marker

3mm 0.75-1mm 2-3 ×l Marker
4-6 ×l 5mm 0.75-1.5mm

Marker

PAGE 0.5-1 ×l Marker 1 loading buffer

