

FS™ Mix	2×	PCR	FS™ Taq
DNA	dNTPs	PCR	
		PCR	
Ī ! ĩ	FS™ Mix	FS™ Taq DNA	
		3'-dA	TA
FS™ Taq DNA		Taq DNA	DNA
	KOD	20s*	

		1		2		
DNA				1: 10000	3	
	DNA	4	PCR	PCR Enhancer	Cat. # P9041	MgCl ₂
	Cat. # P9031					

1			DNA	
2 FS™ DNA			PCR	3' 1
		TA	FS™ Mix	A
3			FS™ Taq DNA	
	1×10 ⁻⁵			

		15-30		3		3
G	C			GC		40-60%
	GC	Tm	55-65		Tm	
				Tm		